

Product datasheet for **MR207688L3V**

Edil3 (NM_001037987) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Edil3 (NM_001037987) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Edil3
Synonyms:	Del-1; Del1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001037987
ORF Size:	1443 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207688).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001037987.3 , NP_001033076.1
RefSeq Size:	5358 bp
RefSeq ORF:	1443 bp
Locus ID:	13612
UniProt ID:	O35474
Cytogenetics:	13 C3



[View online »](#)

Gene Summary:

Promotes adhesion of endothelial cells through interaction with the alpha-v/beta-3 integrin receptor. Inhibits formation of vascular-like structures. May be involved in regulation of vascular morphogenesis of remodeling in embryonic development.[UniProtKB/Swiss-Prot Function]