

Product datasheet for **MR207714L4V**

Foxg1 (NM_001160112) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Foxg1 (NM_001160112) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Foxg1
Synonyms:	2900064B05Rik; BF-1; Bf1; Hfh9; Hfhbf1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001160112
ORF Size:	1443 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207714).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001160112.1 , NP_001153584.1
RefSeq Size:	3986 bp
RefSeq ORF:	1446 bp
Locus ID:	15228
UniProt ID:	Q60987
Cytogenetics:	12 21.78 cM
Gene Summary:	Transcription repression factor which plays an important role in the establishment of the regional subdivision of the developing brain and in the development of the telencephalon. [UniProtKB/Swiss-Prot Function]



[View online »](#)