

## Product datasheet for **MR207723L3V**

### Pde7a (NM\_001122759) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pde7a (NM_001122759) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pde7a
Synonyms:	AU015378; AW047537
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001122759
ORF Size:	1449 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207723).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001122759.1</a> , <a href="#">NP_001116231.1</a>
RefSeq Size:	6298 bp
RefSeq ORF:	1449 bp
Locus ID:	18583
UniProt ID:	<a href="#">P70453</a>
Cytogenetics:	3 5.36 cM
Gene Summary:	Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes. May have a role in muscle signal transduction (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)