

Product datasheet for **MR207730L4V**

Pnliprp2 (NM_011128) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pnliprp2 (NM_011128) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pnliprp2
Synonyms:	PLRP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011128
ORF Size:	1446 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207730).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011128.2 , NP_035258.2
RefSeq Size:	1593 bp
RefSeq ORF:	1449 bp
Locus ID:	18947
UniProt ID:	P17892
Cytogenetics:	19 D2
Gene Summary:	Lipase with broad substrate specificity. Can hydrolyze both phospholipids and galactolipids. Acts preferentially on monoglycerides, phospholipids and galactolipids. Contributes to milk fat hydrolysis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)