

## Product datasheet for **MR207733L4V**

### **Fbxl6 (BC054377) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Fbxl6 (BC054377) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fbxl6
Synonyms:	AU021795; Fbl6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC054377
ORF Size:	1449 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207733).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC054377</a> , <a href="#">AAH54377</a>
RefSeq Size:	1571 bp
RefSeq ORF:	1451 bp
Locus ID:	30840
Cytogenetics:	15 D3
Gene Summary:	Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.[UniProtKB/Swiss-Prot Function]



[View online »](#)