

Product datasheet for **MR207745L3V**

Mecp2 (NM_010788) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mecp2 (NM_010788) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mecp2
Synonyms:	1500041B07Rik; D630021H01Rik; Mbd5; WBP10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010788
ORF Size:	1455 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207745).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010788.3 , NP_034918.1
RefSeq Size:	10576 bp
RefSeq ORF:	1455 bp
Locus ID:	17257
UniProt ID:	Q9Z2D6
Cytogenetics:	X 37.63 cM



[View online »](#)

Gene Summary:

Chromosomal protein that binds to methylated DNA. It can bind specifically to a single methyl-CpG pair. It is not influenced by sequences flanking the methyl-CpGs. Mediates transcriptional repression through interaction with histone deacetylase and the corepressor SIN3. Binds both 5-methylcytosine (5mC) and 5-hydroxymethylcytosine (5hmC)-containing DNA, with a preference for 5-methylcytosine (5mC).[UniProtKB/Swiss-Prot Function]