

Product datasheet for **MR207752L3V**

Osbpl2 (NM_144500) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Osbpl2 (NM_144500) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Osbpl2
Synonyms:	C130070J12Rik; mKIAA0772; ORP-2; Orp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144500
ORF Size:	1455 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207752).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144500.2
RefSeq Size:	2745 bp
RefSeq ORF:	1455 bp
Locus ID:	228983
UniProt ID:	Q8BX94
Cytogenetics:	2 H4



[View online »](#)

Gene Summary:

Intracellular transport protein that binds sterols and phospholipids and mediates lipid transport between intracellular compartments. Increases plasma membrane cholesterol levels and decreases phosphatidylinositol-4,5-bisphosphate levels in the cell membrane. Binds phosphoinositides, such as phosphatidylinositol-4,5-bisphosphate. Exhibits strong binding to phosphatidic acid and weak binding to phosphatidylinositol 3-phosphate. Binds cholesterol, dehydroergosterol, 22(R)-hydroxycholesterol and 25-hydroxycholesterol (in vitro). [UniProtKB/Swiss-Prot Function]