

Product datasheet for **MR207784L4V**

Syt16 (BC056932) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Syt16 (BC056932) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Syt16
Synonyms:	6430703N11, Strep14, syt14r
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC056932
ORF Size:	1458 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207784).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC056932 , AAH56932
RefSeq Size:	2829 bp
RefSeq ORF:	1460 bp
Locus ID:	238266
Cytogenetics:	12 C3
Gene Summary:	May be involved in the trafficking and exocytosis of secretory vesicles in non-neuronal tissues. Is Ca(2+)-independent.[UniProtKB/Swiss-Prot Function]


[View online »](#)