

## Product datasheet for **MR207789L3V**

### Hcls1 (NM\_008225) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hcls1 (NM_008225) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hcls1
Synonyms:	AW213261; HS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008225
ORF Size:	1461 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207789).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008225.1</a> , <a href="#">NP_032251.2</a>
RefSeq Size:	2031 bp
RefSeq ORF:	1461 bp
Locus ID:	15163
UniProt ID:	<a href="#">P49710</a>
Cytogenetics:	16 B3
Gene Summary:	Substrate of the antigen receptor-coupled tyrosine kinase. Plays a role in antigen receptor signaling for both clonal expansion and deletion in lymphoid cells. May also be involved in the regulation of gene expression (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)