

## Product datasheet for **MR207799L3V**

### **Bcl11a (BC010585) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Bcl11a (BC010585) Mouse Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | Bcl11a   |
| Synonyms:                 | 2810047E18Rik; BCL-11A; Ctip1; D930021L15Rik; Evi9; Evi9a; Evi9b; Evi9c; mKIAA1809   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | BC010585   |
| ORF Size:                 | 1461 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR207799).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">BC010585</a> , <a href="#">AAH10585</a>  |
| RefSeq Size:              | 2095 bp  |
| RefSeq ORF:               | 1463 bp  |
| Locus ID:                 | 14025  |
| Cytogenetics:             | 11 A3.2  |
| Gene Summary:             | Functions as a myeloid and B-cell proto-oncogene. May play important roles in leukemogenesis and hematopoiesis. An essential factor in lymphopoiesis, is required for B-cell formation in fetal liver. May function as a modulator of the transcriptional repression activity of ARP1.[UniProtKB/Swiss-Prot Function]  |



[View online »](#)