

Product datasheet for **MR207839L4V**

Krt8 (NM_031170) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Krt8 (NM_031170) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Krt8
Synonyms:	AA960620; AL022697; AU019895; Card2; EndoA; K8; Krt-2.8; Krt2-8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031170
ORF Size:	1473 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207839).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031170.2 , NP_112447.2
RefSeq Size:	1805 bp
RefSeq ORF:	1473 bp
Locus ID:	16691
UniProt ID:	P11679
Cytogenetics:	15 57.2 cM
Gene Summary:	Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.[UniProtKB/Swiss-Prot Function]



[View online »](#)