

## Product datasheet for **MR207869L3V**

### Mns1 (NM\_008613) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mns1 (NM_008613) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mns1
Synonyms:	AW546487
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008613
ORF Size:	1476 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207869).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008613.2</a> , <a href="#">NP_032639.1</a>
RefSeq Size:	1837 bp
RefSeq ORF:	1476 bp
Locus ID:	17427
UniProt ID:	<a href="#">Q61884</a>
Cytogenetics:	9 D
Gene Summary:	May play a role in the control of meiotic division and germ cell differentiation through regulation of pairing and recombination during meiosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)