

## Product datasheet for **MR207893L3V**

### Cnot6l (NM\_144910) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cnot6l (NM_144910) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cnot6l
Synonyms:	4932442K20Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144910
ORF Size:	1653 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207893).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_144910.1</a> , <a href="#">NP_659159.1</a>
RefSeq Size:	8564 bp
RefSeq ORF:	1653 bp
Locus ID:	231464
UniProt ID:	<a href="#">Q8VEG6</a>
Cytogenetics:	5 E3


[View online »](#)

**Gene Summary:**

Poly(A) nuclease with 3'-5' RNase activity. Catalytic component of the CCR4-NOT complex which is one of the major cellular mRNA deadenylases and is linked to various cellular processes including bulk mRNA degradation, miRNA-mediated repression, translational repression during translational initiation and general transcription regulation. Additional complex functions may be a consequence of its influence on mRNA expression. Involved in mRNA decay mediated by the major-protein-coding determinant of instability (mCRD) of the FOS gene in the cytoplasm. Involved in deadenylation-dependent degradation of CDKN1B mRNA. Its mRNA deadenylase activity can be inhibited by TOB1. Mediates cell proliferation and cell survival and prevents cellular senescence (By similarity).[UniProtKB/Swiss-Prot Function]