

Product datasheet for **MR207954L3V**

Oog1 (NM_178657) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Oog1 (NM_178657) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Oog1
Synonyms:	c-1; Oog
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178657
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207954).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178657.4 , NP_848772.3
RefSeq Size:	1840 bp
RefSeq ORF:	1491 bp
Locus ID:	193322
UniProt ID:	E9Q5G7
Cytogenetics:	12 D2
Gene Summary:	May act as a Ras-mediated signaling protein in early embryogenesis.[UniProtKB/Swiss-Prot Function]


[View online »](#)