

## Product datasheet for **MR207961L4V**

### **Pctk1 (BC011069) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pctk1 (BC011069) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pctk1
Synonyms:	Crk5; Pctaire1; Pctk1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC011069
ORF Size:	1488 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207961).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC011069</a> , <a href="#">AAH11069</a>
RefSeq Size:	2815 bp
RefSeq ORF:	1490 bp
Locus ID:	18555
Cytogenetics:	X 16.18 cM



[View online »](#)

**Gene Summary:**

Protein kinase that plays a role in vesicle-mediated transport processes and exocytosis. Can phosphorylate CCNY at 'Ser-336' (in vitro) (By similarity). Plays a role in the regulation of insulin secretion in response to changes in blood glucose levels. Regulates GH1 release by brain neurons. Phosphorylates NSF, and thereby regulates NSF oligomerization. Required for normal spermatogenesis. Regulates neuron differentiation and dendrite development. [UniProtKB/Swiss-Prot Function]