

Product datasheet for **MR207964L3V**

Trim39 (BC031540) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Trim39 (BC031540) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Trim39
Synonyms:	tfp, RBCC-B30.2, mKIAA4179
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC031540
ORF Size:	1488 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207964).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC031540 , AAH31540
RefSeq Size:	2948 bp
RefSeq ORF:	1490 bp
Locus ID:	79263
Cytogenetics:	17 B1
Gene Summary:	E3 ubiquitin-protein ligase. May facilitate apoptosis by inhibiting APC/C-Cdh1-mediated poly-ubiquitination and subsequent proteasome-mediated degradation of the pro-apoptotic protein MOAP1.[UniProtKB/Swiss-Prot Function]



[View online »](#)