

Product datasheet for **MR207977L4V**

Ugp2 (BC025585) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ugp2 (BC025585) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ugp2
Synonyms:	MGC38262
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC025585
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR207977).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC025585 , AAH25585
RefSeq Size:	2069 bp
RefSeq ORF:	1493 bp
Locus ID:	216558
Cytogenetics:	11 13.89 cM
Gene Summary:	Plays a central role as a glucosyl donor in cellular metabolic pathways.[UniProtKB/Swiss-Prot Function]



[View online »](#)