

Product datasheet for **MR207991**

Btbd35f1 (NM_027955) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Btbd35f1
Synonyms:	1700017G21Rik; Gclh; Gmcll1; Gmcl2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >MR207991 representing NM_027955
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC


ATGGGGCTTTTAGTCAGCAGGGTCTTGAGATGCAGGGATTCCAGTCTGTTAGAGCCACAGCCAGAAGCCA
TAGCCGGAGCCAGCTACATTCTGGCAGTCGCAAGCTAAAAAGAAACAGTTTGAGGAGTTGGCAACAAG
TTCTAATGTTTCATGGACCTCAAAACCAGGGAATGTATCCACATCAAGTTCTCACTACATCTACTGGAAG
AGGGTTAAGATCTCATCTAATGATGCTTATCAAACTATTTTTGGACGGGCATGATAGCCACATTAAGAA
TCCGTGCTCTGGGAAAAACATGGTGTTCACAAAGTATTTTTATGTCAGTCAGGCTACTTTGCTAACAT
TCTCAAGGTAAGTCTGGAGAGAATCACACCATGGTGTATAAATCTGATCATTAAAGATGAGGATATTGAT
ACCCGATCTCTGCATTTGTGTTGGTGTCTTGTATACGGATGCGGATTGTCAATAACACCTCTGGAAG
TTCTCAAGTTTTGGCAGCAGCATGCTGCTCGGGTGGATCGAGTAATTCAGCAGTGTGAAGGAATCAT
GAAAGAACTATCAACAGGAACACTGTGTCTCTATTATTTGGCAGCAGAAACCTATAGATTAAAGCT
GTAAAGACGAGATGCTTTGAATGGCTTCTTGCATTTGATGGTACATCCAAGTGTGGCACTTTACAAGG
AAGTAGACTTGAAGTTGATGTATCTCTAGCACTGTCTTCTGACTTACTAGTCATGCAAAAGGAGATTGA
TGTATATACCACTAAAAATATGGATGTTCTTTATCTTAATCCATGCTGGAACGGAACCATGAAACAG
CTTTTACAACAGCAAACTGGCTTTCCACCCACATGGCATATGTTGATAACATCAGTTTTCTTGAAAG
GTGAAGAAGGACTAATATTTCAACCAAGTGTAAAAAGCTGAGATTTAGCAGCATTATCTGTGACTTGAC
TTCCACAACCTATTCTGAACAAGATCGACTAATACCTATGGCATGGTTGTACCCATTTACAACAACAG
TGGTTGACTTTGCTGCAACACAAGAATATGGGGTAATTGGACCACAAGTTATCAATGAACAAGAAGCTTG
AAGAATGCACCATGAGGTGTGATCAATGATCCCAAGGATGGAAGATATACTTGAAGTGGTCAGTTGG
ACGACTTGGCTTTCTTTACGTGTGACCTTTACCAGGCAGTGTGTAATTTAAGGCAACGGTGTGAGAGG
TGTGATGGTTCTGCTTGCACAAACCATATCCGAAATGTCATATTCAGAATAACTTTGGTGTGTTTGATT
CCAACAAAAGAGTAACCTTTCAGAAAGACAAGGTTATAAATCCTCACCTTTGAATATAAGGAGGAGCA
AATTGTAATGAAATGGATAGTGTCTAACCTCCCTATGTGTATATTCTGCAATTTCTTTTGTGTA
AACCTAGGAATGCAGAAAACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/mm9078_g12.zip


Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN:	NM_027955
ORF Size:	1494 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_027955.1 , NM_027955.2 , NM_027955.3 , NP_082231.2
RefSeq Size:	1955 bp
RefSeq ORF:	1497 bp
Locus ID:	71847
Cytogenetics:	X A3.1
MW:	58.2 kDa

Product images:



Circular map for MR207991