

## Product datasheet for **MR208007L4V**

### IL17rb (NM\_019583) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL17rb (NM_019583) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	IL17rb
Synonyms:	Evi27; IL-17ER; IL-17Rh1; IL17br; IL17RH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019583
ORF Size:	1500 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208007).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019583.2</a>
RefSeq Size:	2061 bp
RefSeq ORF:	1500 bp
Locus ID:	50905
UniProt ID:	<a href="#">Q9JIP3</a>
Cytogenetics:	14 18.39 cM
Gene Summary:	Receptor for the proinflammatory cytokines IL17B and IL17E. May play a role in controlling the growth and/or differentiation of hematopoietic cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)