

## Product datasheet for **MR208040L3V**

### **Pwp1 (NM\_133993) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pwp1 (NM_133993) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pwp1
Synonyms:	2310058A11Rik; 2610205J09Rik; AI098170; AL022879; AV083409; AW550690
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133993
ORF Size:	1506 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208040).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133993.1</a> , <a href="#">NP_598754.1</a>
RefSeq Size:	2579 bp
RefSeq ORF:	1506 bp
Locus ID:	103136
UniProt ID:	<a href="#">Q99LL5</a>
Cytogenetics:	10 C1
Gene Summary:	Chromatin-associated factor that regulates transcription (By similarity). Regulates Pol I-mediated rRNA biogenesis and, probably, Pol III-mediated transcription (By similarity). Regulates the epigenetic status of rDNA (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)