

## Product datasheet for **MR208059L3V**

### Sf3a3 (NM\_029157) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Sf3a3 (NM_029157) Mouse Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | Sf3a3  |
| Synonyms:                 | 60kDa; 4930512K19Rik   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_029157  |
| ORF Size:                 | 1506 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR208059).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_029157.3</a> , <a href="#">NP_083433.1</a>  |
| RefSeq Size:              | 1710 bp  |
| RefSeq ORF:               | 1506 bp  |
| Locus ID:                 | 75062  |
| UniProt ID:               | <a href="#">Q9D554</a>   |
| Cytogenetics:             | 4 D2.2   |



[View online »](#)

**Gene Summary:**

Involved in pre-mRNA splicing as a component of the splicing factor SF3A complex that contributes to the assembly of the 17S U2 snRNP, and the subsequent assembly of the pre-spliceosome 'E' complex and the pre-catalytic spliceosome 'A' complex. Involved in pre-mRNA splicing as a component of pre-catalytic spliceosome 'B' complexes.[UniProtKB/Swiss-Prot Function]