

Product datasheet for **MR208126L4V**

Dcaf17 (NM_198005) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dcaf17 (NM_198005) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dcaf17
Synonyms:	2810055O12Rik; 4833418A01Rik; A030004A10Rik; A930009G19Rik; AI448937
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198005
ORF Size:	1518 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208126).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198005.1 , NP_932122.1
RefSeq Size:	1811 bp
RefSeq ORF:	1518 bp
Locus ID:	75763
UniProt ID:	Q3TUL7
Cytogenetics:	2 C2
Gene Summary:	May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex. [UniProtKB/Swiss-Prot Function]



[View online »](#)