

Product datasheet for **MR208160**

Amy2 (BC094924) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Amy2 (BC094924) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Amy2
Synonyms:	1810008N23Rik; Amy-2; Amy2; mAmy2-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR208160 representing BC094924
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGAAGTTCGTTCTGCTGCTTTCCCTCATTGGGTTCTGCTGGGCTCAATATGACCCACATACTTCAGATG
 GGAGGACTGCTATTGTCCACCTGTTTCGAGTGCGCTGGGTTGATATTGCCAAGGAATGTGAGCGATACTT
 AGCTCCTAAGGGATTGAGGGGTGCAGGTCTCTCCACCAATGAAAACGTTGTAGTTTCATAACCCATCA
 AGACCTTGGTGGGAAAGATACCAACCAATCAGCTATAAAATCTGCACAAGGTCTGGAAATGAAGATGAAT
 TCAGAGACATGGTGACAAGGTGCAACAATGTTGGTGTCCGATTATGTGGATGCTGTCATTAACCCAT
 GTGTGGCGCAGGCAATCCTGCAGGAACAAGCAGTACCTGTGGAAGTTACCTCAATCCAAATAACAGGGAA
 TTCCACGAGTTCATACTCTGCTTGGGACTTTAACGATAATAAATGTAATGGAGAAATTGATACTACA
 ATGATGCTTATCAGGTGAGAAATGTCGCTGACTGGCCTTCTGGATCTTGCACTTGAGAAAGATTATGT
 TCGTACCAAGGTGGCTGACTATATGAACCATCTCATTGACATTGGAGTAGCAGGGTTCAGACTTGATGCT
 GCTAAGCACATGTGGCCTGGAGACATAAAGGCAGTTTTGGACAACTGCATAATCTCAATACAAAATGGT
 TCTCCCAAGGAAGCAGACCTTTCAATTTCCAAGAGGTCATTGATCTGGGTGGTGAGGCAATTAAAGGTAG
 TGAGTACTTTGGAATGGCCGTGTGACAGAATCAAGTATGGTGCAAACTTGGCACAGTTATCCGCAAG
 TGGAATGGCGAGAAGATGTCCTATTTAAGAAGTGGGAGAGGTTGGGTTTGGTGCCTTCTGACAGAG
 CCCTTGTGTTTGTGGACAACCATGATAATCAGCGAGGACATGGTCTGGAGGATCATCCATCCTGACATT
 CTGGGATGCTAGAATGTATAAAATGGCTGTGCGATTTATGTTGGCTCATCCTTATGGATTACAAGAGTA
 ATGTCAAGTTACCGTTGGAATAGAAATTTCCAGAATGGAAGATCAGAATGACTGGATTGGACACCCA
 ATAACAATGGAGTAACAAAAGAAGTGACCATTAATGCAGACACTACTTGTGGCAATGACTGGGTCTGTGA
 ACACAGATGGCGTCAAATCAGGAACATGGTTGCCCTTCAGGAATGTGGTCAATGGTCAGCCTTTTCAAAC
 TGGTGGGATAATAACAGCAACCAAGTAGCTTTAGCAGAGGAAACAGAGGATTCAATGTCTTTAAACAATG
 ATGACTGGGCTTTGTGAGCCACTTTACAGACTGGTCTTCTGCTGGCACATACTGTGATGTCATCTCTGG
 AGATAAGGTGATGGCAATTGCACTGGACTTAGAGTGAATGTTGGCAGTGATGGCAAAGCTCACTTTTCC
 ATTAGTAACCTGCTGAGGACCCATTTATTGCAATCCATGCTGACTCAAAATTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR208160 representing BC094924
 Red=Cloning site Green=Tags(s)


MKFLVLLSLIGFCWAQYDPHTSDGRTAIVHLFEWRWVDIAKECERYLAPKGFGGVQVSPNENNVVHNPS
 RPWWERYQPISYKICTRSGNEDEFDMVTRCANNVGVRIYVDVINHMCGAGNPAGTSSTCGSYLNPNNRE
 FPAVPYSAWDFNDKNCNGEIDNYNDAYQVRNCRLTGLLDLALEKDYVRTKVADYMNHLIDIGVAGFRLDA
 AKHMWPGDIKAVLDKHLNLTWFSQGSRPFIQFEVIDLGGEAIKGYEYFGNGRVTEFKYGAKLGTVIRK
 WNGEKMSYLKNWEGEWGLVPSDRALVFVDNHDNRGHGAGGSSILTFWDARMYKMAVGFMALHPYGFTRV
 MSSYRWRNRFQNGKDQNDWIGPPNNNGVTKEVTINADTTGNDWVCEHRWRQIRNMVAFRNVVNGQPF SN
 WWDNNSNQVAFSRGNRGFIVFNDDWALSATLQTGLPAGTYCDVISGDKVDGNGCTGLRVNVGSDGKAHFS
 ISNSAEDPFIAIHADSKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: BC094924

ORF Size: 1524 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC094924 , AAH94924
RefSeq Size:	1589 bp
RefSeq ORF:	1526 bp
Locus ID:	11723
MW:	58.2 kDa

Product images:



Circular map for MR208160