

## Product datasheet for **MR208161L3V**

### **Edc3 (NM\_153799) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Edc3 (NM_153799) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Edc3
Synonyms:	AA517853; Lsm16; Yjdc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153799
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208161).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153799.1</a> , <a href="#">NP_722494.1</a>
RefSeq Size:	3343 bp
RefSeq ORF:	1527 bp
Locus ID:	353190
UniProt ID:	<a href="#">Q8K2D3</a>
Cytogenetics:	9 B
Gene Summary:	Binds single-stranded RNA. Involved in the process of mRNA degradation and in the positive regulation of mRNA decapping (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)