

Product datasheet for **MR208163L3V**

Frs2 (NM_177798) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Frs2 (NM_177798) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Frs2
Synonyms:	4732458E18; C330018A15Rik; Frs2alpha; SNT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177798
ORF Size:	1524 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208163).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177798.3 , NP_808466.1
RefSeq Size:	5701 bp
RefSeq ORF:	1527 bp
Locus ID:	327826
UniProt ID:	Q8C180
Cytogenetics:	10 D2



[View online »](#)

Gene Summary:

Adapter protein that links activated FGR and NGF receptors to downstream signaling pathways. Plays an important role in the activation of MAP kinases and in the phosphorylation of PIK3R1, the regulatory subunit of phosphatidylinositol 3-kinase, in response to ligand-mediated activation of FGFR1. Modulates signaling via SHC1 by competing for a common binding site on NTRK1.[UniProtKB/Swiss-Prot Function]