

Product datasheet for **MR208164L4V**

Tle6 (BC048832) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tle6 (BC048832) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tle6
Synonyms:	Grg6, 1810057E06Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC048832
ORF Size:	1524 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208164).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC048832 , AAH48832
RefSeq Size:	1897 bp
RefSeq ORF:	1526 bp
Locus ID:	114606
Cytogenetics:	10 C1
Gene Summary:	As a member of the subcortical maternal complex (SCMC), plays an essential role for zygotes to progress beyond the first embryonic cell divisions.[UniProtKB/Swiss-Prot Function]



[View online »](#)