

Product datasheet for **MR208186L4V**

Cttn (BC011434) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cttn (BC011434) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cttn
Synonyms:	1110020L01Rik; Ems1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC011434
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208186).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC011434 , AAH11434
RefSeq Size:	2916 bp
RefSeq ORF:	1529 bp
Locus ID:	13043
Cytogenetics:	7 F5



[View online »](#)

Gene Summary:

Contributes to the organization of the actin cytoskeleton and cell shape (PubMed:17403031). Plays a role in the formation of lamellipodia and in cell migration (By similarity). Plays a role in the regulation of neuron morphology, axon growth and formation of neuronal growth cones (By similarity). Through its interaction with CTTNBP2, involved in the regulation of neuronal spine density (PubMed:22262902). Plays a role in the invasiveness of cancer cells, and the formation of metastases (By similarity). Plays a role in focal adhesion assembly and turnover (By similarity). In complex with ABL1 and MYLK regulates cortical actin-based cytoskeletal rearrangement critical to sphingosine 1-phosphate (S1P)-mediated endothelial cell (EC) barrier enhancement (By similarity). Plays a role in intracellular protein transport and endocytosis, and in modulating the levels of potassium channels present at the cell membrane (PubMed:17959782). Plays a role in receptor-mediated endocytosis via clathrin-coated pits (By similarity). Required for stabilization of KCNH1 channels at the cell membrane (By similarity).[UniProtKB/Swiss-Prot Function]