

Product datasheet for **MR208196L3V**

C230094A16Rik (BC024726) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C230094A16Rik
Synonyms:	MGC38936, RP23-176O5.5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC024726
ORF Size:	1530 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR208196).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [BC024726](#), [AAH24726](#)

RefSeq Size: 2444 bp

RefSeq ORF: 1532 bp

Locus ID: 237711

Cytogenetics: 11 A3.3

Gene Summary: May modify the assembly dynamics of microtubules, such that microtubules are slightly longer, but more dynamic.[UniProtKB/Swiss-Prot Function]



