

Product datasheet for **MR208203L3V**

Etv5 (NM_023794) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Etv5 (NM_023794) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Etv5
Synonyms:	1110005E01Rik; 8430401F14Rik; ERM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023794
ORF Size:	1533 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208203).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023794.2 , NP_076283.2
RefSeq Size:	3847 bp
RefSeq ORF:	1533 bp
Locus ID:	104156
UniProt ID:	Q9CXC9
Cytogenetics:	16 B1
Gene Summary:	Binds to DNA sequences containing the consensus nucleotide core sequence 5'-GGAA.-3'. [UniProtKB/Swiss-Prot Function]



[View online »](#)