

Product datasheet for **MR208207L4V**

Aldh7a1 (NM_001127338) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aldh7a1 (NM_001127338) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aldh7a1
Synonyms:	Atq1; D18Wsu181e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001127338
ORF Size:	1536 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208207).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001127338.1 , NP_001120810.1
RefSeq Size:	3105 bp
RefSeq ORF:	1536 bp
Locus ID:	110695
UniProt ID:	Q9DBF1
Cytogenetics:	18 30.55 cM



[View online »](#)

Gene Summary:

Multifunctional enzyme mediating important protective effects. Metabolizes betaine aldehyde to betaine, an important cellular osmolyte and methyl donor. Protects cells from oxidative stress by metabolizing a number of lipid peroxidation-derived aldehydes. Involved in lysine catabolism (By similarity).[UniProtKB/Swiss-Prot Function]