

Product datasheet for **MR208220L4V**

Oasl1 (NM_145209) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Oasl1 (NM_145209) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Oasl1
Synonyms:	7530414C13Rik; oasl9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145209
ORF Size:	1536 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208220).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145209.2 , NP_660210.1
RefSeq Size:	2124 bp
RefSeq ORF:	1536 bp
Locus ID:	231655
UniProt ID:	Q8VI94
Cytogenetics:	5 F
Gene Summary:	Does not have 2'-5'-OAS activity, but can bind double-stranded RNA. Displays antiviral activity via an alternative antiviral pathway independent of RNase L.[UniProtKB/Swiss-Prot Function]



[View online »](#)