

Product datasheet for **MR208252L4V**

Gtf3c5 (NM_148928) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gtf3c5 (NM_148928) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gtf3c5
Synonyms:	2700084A09Rik; TFiiC2-63; TFIIIC63
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_148928
ORF Size:	1545 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208252).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_148928.1
RefSeq Size:	2805 bp
RefSeq ORF:	1545 bp
Locus ID:	70239
UniProt ID:	Q8R2T8
Cytogenetics:	2 A3
Gene Summary:	Involved in RNA polymerase III-mediated transcription. Integral, tightly associated component of the DNA-binding TFIIIC2 subcomplex that directly binds tRNA and virus-associated RNA promoters (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)