

Product datasheet for **MR208255L4V**

Baiap2l1 (NM_025833) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Baiap2l1 (NM_025833) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Baiap2l1
Synonyms:	1300006M19Rik; AI585895; IRTKS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025833
ORF Size:	1542 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208255).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025833.3 , NP_080109.1
RefSeq Size:	3248 bp
RefSeq ORF:	1545 bp
Locus ID:	66898
UniProt ID:	Q9DBJ3
Cytogenetics:	5 G2
Gene Summary:	May function as adapter protein. Involved in the formation of clusters of actin bundles. Plays a role in the reorganization of the actin cytoskeleton in response to bacterial infection (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)