

## Product datasheet for MR208256L3

### Tmem117 (NM\_178789) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem117 (NM_178789) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Tmem117
Synonyms:	B930062P21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208256).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_178789
ORF Size:	1542 bp



[View online »](#)

