

Product datasheet for **MR208284**

2010005J08Rik (BC066994) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	2010005J08Rik (BC066994) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	2010005J08Rik
Synonyms:	2010005J08Rik; 9130001I21Rik; AW060359; mKIAA1507; Urg4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR208284 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAAGGAGATGACTGCGGGTTCATTATGGAGATGGTACAAACGAGGCTCAGGACAATGACTTCCCCA
 CAGTGGAGAGAAGCAGGCTTCAGGAGATGCTGTCCCTGTTGGGACTGGAGACATACCAGGCACAGAAGCT
 CACCCTCCAGGACTCCCTGCAGATCAGTTTTGACAGCATGAAGAACTGGGCTCCTCAGGTTCCCAAGGAC
 CTTCCCTGGCATTTCCTCAGGAAGCTGCAGGCACTCAACGCCGAAGCCAGGAACACCACCATGGTCTGG
 ACGTGCCTCCAGAGGCTGGCTGCAGAGAAGGAGAGCCAGGCAGAGGAGGAGATGATGTACTGGGACCC
 AGCGGAGGACGTCGCTGCTGACATCTACTCTTTCTCAGAGCTGCCATGCCTGATACGCCTGTGAACCCC
 CTGGATCTTCTGCGCCCTCTGCTGTCTTCAGACACCTTTCTGCAGCAGGAAGTCTGCTGAAGATGT
 CCCTCTGTCAGTTTGCACTCCACTTGTGTTGCCCGATTCTGAAAACCACTACCACACCTTTCTGCTCTG
 GGCCTGCGGGGTGTTGTGCGGACGTGGTGGTTCCAGCCCTCCGGAGCCTGGGGAGTTTCCGAGAAGAC
 AGTGTGGTCTGTCCAGGGTGCCACGTTTGCCTTCGTGCGCATGGATGTCAGCAGCAACTCCAAGTCCC
 AGCTGCTCAATGCCGTCTCAGCCCGGCCCGCCAGTGGGACTGCTTCTGGCATCGGGACCTCAATTT
 AGGCACTAATCCTCGCGAGATCTCGGATGGGTTAGTAGAGATCTCCTGGTTTTTCCCCAGTGGTAAGGAG
 GACCTGGATGTTTTCCGGAGCCCATGGCTTTTGAACCTTGAGAGGTGACATTGGGTCTCACTGGCTGC
 AGTTCAAGCTCTTGACAGAAGTCTCCTCTGCTGTGTTATTTTACTGACAACATCAGCAAGAAGGAATA
 CAAACTGCTGTACTCCATGAAGGAGTCAACCACAAAGTACTACTTCATCCTGAGTCCCTACCGAGGGAAA
 CGGAACACGAACCTGCGCTTCTCAACAGGCTGATCCCAGTGTGAAGATAGACCCTCTCACGTGCTGG
 TGAAGTTCAGCAGCACGGACAGCGAGGGCTTCGTGAGGCGCATCGAGCCATCATGTGACTGTGGCTCG
 CTCTCCCTGCAGGAGGGTCTGTGGAGGACATGGCACACGCGGCACGTAAGCTCGGCCTCAGAGTTGAT
 GAGGACTTTGAGGAGTGTGAGAAGGCAAAGGACCGGATGGAGAAAATCCTCAGGAAGATCAAGGATGTGG
 ACGCCTACAAGAGAGACGAGCTGCGGCTGCAGGGCGAGCCTGCGAGGAGGGCCGCCAAGCTGAAAGAGA
 GTTCTGCCAGCTCCAGTGGGTCCCGAGCCCCAGAGAAGCACCGGGCCGAGCTGAGACGGCGCTCCTT
 GAGCTCCGGGTGCAGCAGAATGGCCAGGAGCCACCTCAGGGTCCAGGAGTTCATCTGGGAATAAGCA
 GCCCTCCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR208284 protein sequence
 Red=Cloning site Green=Tags(s)

MEGDDCGFHYDGTNEAQDNDPPTVERSRLQEMLSLGLLETYQAQKLTLDLSLQISFDSMKNWAPQVPKD
 LPWHFLRKLQALNAEARNTTMVLDVPEAWPAEKESQAEEEMMYWPAEDVAADIYSFSELPMPDTPVNP
 LDLLCALLLSSDFFLQQEVLLKMSLCQFALPLVLPDSENHYHTFLLWALRGVVRTWWFQPLRSLGSRFED
 SVVLSRVPTFAFVRMDVSSNSKSQLLNVLSPARRQWDCFWRDLNLTGTPREISDGLVEISWFFPSGKE
 DLDFVPEPMAFLNLRGDI GSHWLQFKLLTEVSSAVFILTDNISKEYKLLYSMKESTTKYFILSPYRGK
 RNTNLRFLNRLIPVLKIDHSHLVKVSSTDSEGFVRRIRAIMCSVARSPCRRVSVEDMAAARKLGLRVD
 EDFEECQKAKDRMEKILRKIKVDAYKRDELRLQGEPARRAAQAEREFQQLQWVPEPPEKHRAELRRRVL
 ELRVQNGQEPTSGVQEFILGISSPSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: BC066994

ORF Size: 1551 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC066994](#), [AAH66994](#)

RefSeq Size: 3342 bp

RefSeq ORF: 1553 bp


Locus ID: 72046

Cytogenetics: 11 A1

MW: 59.7 kDa

Gene Summary: May be involved in cell cycle progression through the regulation of cyclin D1 expression. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR208284