

Product datasheet for **MR208284L3V**

2010005J08Rik (BC066994) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	2010005J08Rik (BC066994) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	2010005J08Rik
Synonyms:	2010005J08Rik; 9130001I21Rik; AW060359; mKIAA1507; Urg4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC066994
ORF Size:	1551 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208284).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC066994 , AAH66994
RefSeq Size:	3342 bp
RefSeq ORF:	1553 bp
Locus ID:	72046
Cytogenetics:	11 A1
Gene Summary:	May be involved in cell cycle progression through the regulation of cyclin D1 expression. [UniProtKB/Swiss-Prot Function]



[View online »](#)