

## Product datasheet for **MR208304L4V**

### Camk2g (NM\_001039138) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Camk2g (NM_001039138) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Camk2g
Synonyms:	Camkg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001039138
ORF Size:	1557 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208304).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001039138.1</a> , <a href="#">NP_001034227.1</a>
RefSeq Size:	3693 bp
RefSeq ORF:	1557 bp
Locus ID:	12325
UniProt ID:	<a href="#">Q923T9</a>
Cytogenetics:	14 A3



[View online »](#)

**Gene Summary:**

Calcium/calmodulin-dependent protein kinase that functions autonomously after Ca(2+)/calmodulin-binding and autophosphorylation, and is involved in sarcoplasmic reticulum Ca(2+) transport in skeletal muscle and may function in dendritic spine and synapse formation and neuronal plasticity. In slow-twitch muscles, is involved in regulation of sarcoplasmic reticulum (SR) Ca(2+) transport and in fast-twitch muscle participates in the control of Ca(2+) release from the SR through phosphorylation of the ryanodine receptor-coupling factor triadin. In neurons, may participate in the promotion of dendritic spine and synapse formation and maintenance of synaptic plasticity which enables long-term potentiation (LTP) and hippocampus-dependent learning (By similarity).[UniProtKB/Swiss-Prot Function]