

Product datasheet for **MR208328L3V**

Ptbp3 (NM_144904) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Ptbp3
Synonyms:	5830471K22Rik; AA407443; AI462022; AW107884; C86549; Rod1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PSI00092)
Tag:	Myc-DDK
ACCN:	NM_144904
ORF Size:	1560 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR208328).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_144904.2
RefSeq Size:	6880 bp
RefSeq ORF:	1656 bp
Locus ID:	230257
UniProt ID:	Q8BHD7
Cytogenetics:	4 B3



Gene Summary:

RNA-binding protein that mediates pre-mRNA alternative splicing regulation. Plays a role in the regulation of cell proliferation, differentiation and migration. Positive regulator of EPO-dependent erythropoiesis. Participates in cell differentiation regulation by repressing tissue-specific exons. Promotes Fas exon 6 skipping. Binds RNA, preferentially to both poly(G) and poly(U) (By similarity).[UniProtKB/Swiss-Prot Function]