

## Product datasheet for MR208335L3

### Arfgap2 (NM\_023854) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Arfgap2 (NM_023854) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Arfgap2
Synonyms:	2310032E02Rik; Zfp28; Zfp289
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208335).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_023854
ORF Size:	1563 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_023854.1</a> , <a href="#">NP_076343.2</a>
<b>RefSeq Size:</b>	3367 bp
<b>RefSeq ORF:</b>	1563 bp
<b>Locus ID:</b>	77038
<b>UniProt ID:</b>	<a href="#">Q99K28</a>
<b>Cytogenetics:</b>	2 E1
<b>Gene Summary:</b>	This gene encodes a zinc-finger-containing GTPase-activating protein for ADP ribosylation factor 1 (ARF1), a small GTPase that plays a role in coatamer-mediated vesicular trafficking. This gene product stimulates the hydrolysis of ARF1-bound GTP, which may lead to dissociation of coatamer from Golgi-derived membranes to allow fusion with target membranes. It may regulate the retrograde transport from the Golgi complex to the endoplasmic reticulum. Expression of this gene has been shown to be controlled by inhibitor of DNA binding 1 (Id1). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. A pseudogene of this gene was identified on chromosome 6. [provided by RefSeq, Oct 2009]

Product images:



Circular map for MR208335L3