

Product datasheet for **MR208335L3V**

Arfgap2 (NM_023854) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arfgap2 (NM_023854) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arfgap2
Synonyms:	2310032E02Rik; Zfp28; Zfp289
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023854
ORF Size:	1563 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208335).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023854.1 , NP_076343.2
RefSeq Size:	3367 bp
RefSeq ORF:	1563 bp
Locus ID:	77038
UniProt ID:	Q99K28
Cytogenetics:	2 E1



[View online »](#)

Gene Summary:

This gene encodes a zinc-finger-containing GTPase-activating protein for ADP ribosylation factor 1 (ARF1), a small GTPase that plays a role in coatamer-mediated vesicular trafficking. This gene product stimulates the hydrolysis of ARF1-bound GTP, which may lead to dissociation of coatamer from Golgi-derived membranes to allow fusion with target membranes. It may regulate the retrograde transport from the Golgi complex to the endoplasmic reticulum. Expression of this gene has been shown to be controlled by inhibitor of DNA binding 1 (Id1). Alternatively spliced transcript variants encoding different isoforms have been found for this gene. A pseudogene of this gene was identified on chromosome 6. [provided by RefSeq, Oct 2009]