

Product datasheet for **MR208352L4V**

Slc18a1 (NM_153054) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc18a1 (NM_153054) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc18a1
Synonyms:	4832416I10Rik; Vat1; Vmat1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153054
ORF Size:	1566 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208352).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153054.2
RefSeq Size:	3118 bp
RefSeq ORF:	1566 bp
Locus ID:	110877
UniProt ID:	Q8R090
Cytogenetics:	8 33.88 cM
Gene Summary:	Involved in the transport of biogenic monoamines, such as serotonin, from the cytoplasm into the secretory vesicles of neuroendocrine and endocrine cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)