

Product datasheet for **MR208385L4V**

Tmem63b (BC057136) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem63b (BC057136) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem63b
Synonyms:	BC026370
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC057136
ORF Size:	1569 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208385).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC057136 , AAH57136
RefSeq Size:	1702 bp
RefSeq ORF:	1571 bp
Locus ID:	224807
Cytogenetics:	17 B3
Gene Summary:	Acts as an osmosensitive calcium-permeable cation channel.[UniProtKB/Swiss-Prot Function]


[View online »](#)