

## Product datasheet for **MR208404L3V**

### Arfgap3 (NM\_025445) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Arfgap3 (NM_025445) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arfgap3
Synonyms:	0610009H19Rik; 1810004P07Rik; 1810035F16Rik; 9130416J18Rik; Arfgap1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025445
ORF Size:	1575 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208404).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025445.2</a> , <a href="#">NP_079721.2</a>
RefSeq Size:	2618 bp
RefSeq ORF:	1575 bp
Locus ID:	66251
UniProt ID:	<a href="#">Q9D8S3</a>
Cytogenetics:	15 E1
Gene Summary:	GTPase-activating protein (GAP) for ADP ribosylation factor 1 (ARF1). Hydrolysis of ARF1-bound GTP may lead to dissociation of coatomer from Golgi-derived membranes to allow fusion with target membranes (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)