

## Product datasheet for **MR208408L4V**

### Arhgef3 (BC005517) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgef3 (BC005517) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgef3
Synonyms:	1200004I24Rik; 9830169H03Rik; C76747; Xpln
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC005517
ORF Size:	1572 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208408).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC005517.1</a>
RefSeq Size:	3628 bp
RefSeq ORF:	1574 bp
Locus ID:	71704
Cytogenetics:	14 16.09 cM
Gene Summary:	Acts as guanine nucleotide exchange factor (GEF) for RhoA and RhoB GTPases. [UniProtKB/Swiss-Prot Function]



[View online »](#)