

Product datasheet for **MR208413L4V**

Gmcl1 (NM_011818) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gmcl1 (NM_011818) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gmcl1
Synonyms:	2810049L19Rik; Btbd13; DIP; Gcl; Gcl-1; glc-1; mgcl-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011818
ORF Size:	1572 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208413).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011818.3 , NP_035948.3
RefSeq Size:	2967 bp
RefSeq ORF:	1575 bp
Locus ID:	23885
UniProt ID:	Q920G9
Cytogenetics:	6 37.75 cM
Gene Summary:	Possible function in spermatogenesis. Enhances the degradation of MDM2 and increases the amount of p53 probably by modulating the nucleocytoplasmic transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)