

Product datasheet for **MR208431L3V**

Slc32a1 (NM_009508) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc32a1 (NM_009508) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc32a1
Synonyms:	R75019; VGAT; Viaat
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009508
ORF Size:	1578 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208431).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009508.2 , NP_033534.2
RefSeq Size:	2805 bp
RefSeq ORF:	1578 bp
Locus ID:	22348
UniProt ID:	O35633
Cytogenetics:	2 H1
Gene Summary:	Involved in the uptake of GABA and glycine into the synaptic vesicles.[UniProtKB/Swiss-Prot Function]



[View online »](#)