

Product datasheet for **MR208451L4V**

Oxsr1 (NM_133985) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Oxsr1 (NM_133985) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Oxsr1
Synonyms:	2210022N24Rik; 2810422B09Rik; AI462649; AW209236; mKIAA1101; Osr1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133985
ORF Size:	1581 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208451).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133985.1 , NP_598746.2
RefSeq Size:	1921 bp
RefSeq ORF:	1584 bp
Locus ID:	108737
UniProt ID:	Q6P9R2
Cytogenetics:	9 F3
Gene Summary:	Phosphorylates RELL1, RELL2, RELT and PAK1. Phosphorylates PLSCR1 in the presence of RELT.[UniProtKB/Swiss-Prot Function]



[View online »](#)