

Product datasheet for **MR208467L4V**

Ugt2a2 (NM_001024148) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ugt2a2 (NM_001024148) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ugt2a2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001024148
ORF Size:	1584 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208467).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024148.1 , NP_001019319.1
RefSeq Size:	2595 bp
RefSeq ORF:	1587 bp
Locus ID:	552899
UniProt ID:	Q6PDD0
Cytogenetics:	5 E1
Gene Summary:	UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds.[UniProtKB/Swiss-Prot Function]



[View online »](#)