

Product datasheet for MR208489L3

Traf6 (NM_009424) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Traf6 (NM_009424) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Traf6
Synonyms:	2310003F17Rik; AI851288; C630032O20Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208489).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_009424
ORF Size:	1593 bp




[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_009424.2 , NP_033450.2
RefSeq Size:	5536 bp
RefSeq ORF:	1593 bp
Locus ID:	22034
UniProt ID:	P70196
Cytogenetics:	2 E2
Gene Summary:	This gene encodes a member of the TNF receptor associated factor (TRAF) family of adaptor proteins that mediate signaling events from members of the TNF receptor and Toll/IL-1 receptor families to activate transcription factors such as NF-kappa-B and AP-1. The product of this gene is essential for perinatal and postnatal survival. Mice deficient in this protein exhibit osteopetrosis and defective in development of epidermal appendices, normal B cell differentiation, lymph node organogenesis, interleukin-1 signaling, lipopolysaccharide signaling and neural tube closure. This protein possesses ubiquitin ligase activity. Alternate splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2014]

Product images:



Circular map for MR208489L3



Double digestion of MR208489L3 using SgfI and MluI