

Product datasheet for **MR208499L3V**

Csgalnact1 (NM_172753) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Csgalnact1 (NM_172753) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Csgalnact1
Synonyms:	4732435N03Rik; CSGalNAct-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172753
ORF Size:	1590 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208499).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172753.5 , NP_766341.4
RefSeq Size:	3943 bp
RefSeq ORF:	1593 bp
Locus ID:	234356
UniProt ID:	Q8BJQ9
Cytogenetics:	8 B3.3



[View online »](#)

Gene Summary:

Transfers 1,4-N-acetylgalactosamine (GalNAc) from UDP-GalNAc to the non-reducing end of glucuronic acid (GlcUA). Required for addition of the first GalNAc to the core tetrasaccharide linker and for elongation of chondroitin chains. Important role in chondroitin chain biosynthesis in cartilage formation, and subsequent endochondral ossification (PubMed:17145758, PubMed:21148564). Moreover, is involved in the metabolism of aggrecan (PubMed:21148564).[UniProtKB/Swiss-Prot Function]